



# LABORATORY

# MICROSCOPE



## SW/XSZ-107 207 SERIES BIOLOGICAL MICROSCOPE

OUTFIT	MODEL						
	SW/XSZ -107BN	SW/XSZ -N107	SW /XSZ -N107CCD	SW/XSZ -207	SW/XSZ -N207	SW/XSZ -N207CCD	
Viewing Head	Monocular Head, at 45° Rotatable						
	Sliding Binocular Head, at 45° Rotatable	●	●	●			
	Compensation Free Binocular Head, at 30° Rotatable				●	●	●
	Trinocular Head for Photography						
	Trinocular Head for Video Camera (C Mount)			●			
Objective	Achromatic DIN 4X, 10X, 40X(S), 100X(S, Oil)	●	●	●	●	●	
Eyepiece	Wide Field Eyepiece WF10X, P16X	●	●	●	●	●	
Stage	Double Layers Mechanical Stage 140X140mm/75X75mm		●	●	●	●	
	Single Layers Mechanical Stage 140X140mm/75X75mm	●					
Condenser	Abbe NA1.25 with Iris Diaphragm & Filter	●	●	●	●	●	
Illumination	Built-in Illumination, Halogen Lamp 110V or 220V/6V/20W	●	●	●	●	●	
	Plan-concave Mirror	●	●	●	●	●	



## SW/XSZ-20 SERIES BIOLOGICAL MICROSCOPE SW/XSZ-2001/2002

OUTFIT	MODEL			
	2001	2001A	2002	2002A
Achromatic Objectives	4X. 10X. 40X ( S ).100X( S, Oil)			
Wide Field Eyepiece	WF10X. WF16X			
Sliding Binocular Head	Inclined at 45° , 360° Rotatable Interpupillary Distance 55~75mm			
Nosepiece	Quadruple			
	Quintuple			
Double layer Mechanical stage	Size 140x140mm, Moving Range 75x45mm			
Focusing	Coaxial Coarse and Fine Adjustment, Focusing Range 30mm, Focusing Interval 0.002mm			
Illumination, condenser NA=1.25	6V/20W Halogen Lamp, Adjustable Brightness			
	Kohler Illumination			
Power	AC Input 220V or 110V, DC Output 6V			



## MONOCULAR BIOLOGICAL MICROSCOPE


OUTFIT	MODEL				
	SW /XSZ -2009	SW/XSP -101	SW/XSP -102	SW/XSP -103	SW/SME -F
Achromatic Objectives	4X. 10X. 40X ( S )				
	100X( S, Oil)				
Eyepiece	P10X				
	H5X. H10X. H16X				
	WF10x				
Nosepiece	Triple				
	Quadruple				
Stage	Plain Stage with Slide Clips				
	Double Layer Mechanical Stage				
Focusing	Separate Coarse				
	Coaxial Coarse				
Condenser	Single Lens with Disc Diaphragm				
	Abbe NA1.25 Condenser				
Illumination	Plan-concave Mirror				
	Incandescent Lamp 220V/20W or 110V/20W				
	Halogen Lamp 6V/20W				
Power	AC Input 220V or 110V				

# LABORATORY PRODUCTS



## SW7101, 7102, 7105 SLIDE

Ground edges  
Ungrounded edges  
Frosted one end, ground edges




## SW7102 COVER SLIP

Size: 18 x 18mm, 20 x 20mm, 22 x 22mm, 24 x 24mm,  
24 x 32mm, 24 x 40mm, 24 x 50mm, 24 x 60mm, etc.



## SW715 URIC TEST STRIPS SW716 BLOOD TEST STRIPS



## SW752 RAPID DIAGNOSTIC TESTS

HCG Midstream, HBsAg, HCV, HIV, Morphine, etc.



## BLOOD LANCET SERIES SW7351 FLAKE TYPE SW7352 TWIST TYPE SW7353 CAPPED TYPE



## SW7355 SAFETY TYPE SW7356 AUTO-SAFETY TYPE



## SW738 VACUUM TUBE

Glass tube or plastic tube available



## SW802 ULTRASONIC GEL

250ml, 2.5L, 5L etc.

# LABORATORY PRODUCTS



**SW730 APPLICATOR**  
Length: 3", 6", 8", etc  
Material: wooden, plastic  
With or without cotton tip



**SW731 PIPETTOR**  
Fixed volume  
Adjustable volume



**SL3 SERIES CENTRIFUGETUBE**



**SL4 SERIES PIPETTE TIPS**



**SL5 SERIES CULTURE PLATES & MICROPLATES**



**SL6 SERIES PASTEUR PIPETTES**



**HIGH SPEED MICROHEMATOCRIT CENTRIFUGE**



**SW732 TABLE TYPE CENTRIFUGE**

	80-1	80-2	90-1
MAX. SPEED	4000r/min	4000r/min	5000r/min
MAX. RCF.	1790X g	2130X g	3130Xg
POWER	25w	40w	40w

	SH120-1	SH120-2
MAX.SPEED	12000r/min	12000r/min
MAX.RCF.	15300Xg	15300Xg
POWER OUTPUT	100W	350W
CAPACITY	24PCS CAPILLARY ø1.5mm	24PCS CAPILLARY ø3mm
TIME RANGE	0~10min	0~99min
N.WT.	11.5KGS	11.5KGS
REGULATION	STEP	STEPLESS



**SW820 MEDICAL PACKING**

- Self-sealing flat pouch
- Heat-sealing flat pouch
- Heat-sealing gusseted pouch
- Heat-sealing flat reel
- Heat-sealing gusseted reel
- Self-sealing connecting pouches
- Heat-sealing connecting pouches



**SW804 DISPOSABLE SURGICAL RAZOR SINGLE OR DOUBLE LAYER**



**SW805 DIAGNOSTIC PENLIGHT**




**SW808 WOODEN TONGUE DEPRESSOR**

- Size: Adult, Child
- Material: wooden, plastic
- Individual sterile packing available



**SW809 SURGICAL BRUSH**



**SW810 DISPOSABLE SURGICAL BRUSH**



**SW812 TEST TUBE BRUSH**



**SW815 DISPOSABLE FORCEPS**